

Since 1935
CHIRANA - European quality




Made in Slovakia





-  Plastic parts, electronic CHIRANA, micromotors CHIRANA, metallic parts
-  Amalgam separator, raw material, machine equipment,
-  Micromotors, scalers, upholstery for patient chair, raw material, machine equipment
-  OP lights, curing lamp, valves, syringes, suction systems
-  Driving units for patient chairs
-  Machine equipment
-  Aluminium casting parts
-  Scalers
-  Scalers

CHEESE EASY



Spittoon mounted on patients chair or to floor



Dentist tables

Dentist tables 5N on pantographic arm



Cart dentist tables 5N



Dentist buttons control panel



- Spray and power set up
- Light for instrument on/off
- Chip blower function on/off
- Reverse of MM on/off
- Giromatic function on/off
- Mode control (by foot control or manual from keyboard)
- Disinfection programm - activation
- Glass filling
- Spittoon rinsing
- Stop function
- OP light switch on/off
- Chair control

System of instruments tubes leading

Whip system



Hanging system



Additional tray for dentist table

fixed beside dentist table



Under dentist table

- 1x silicone tray
- 1x stainless tray
- 2x silicone tray
- 2x stainless tray



Spittoon block

Suction systems:

- Ejector saliva
- Ejector suction
- CATTANI
- DURR MSBV
- DURR CS1
- DURR CAS1

Hygiene:

- Clean water bottle
- Disinfection system
- Clean water bottle + disinfection system



Assistant tables


2N – two places



4N – four places



Control panel on assistant table



- Glass filling function
- Spittoon rinsing
- Stop function
- OP light switch on/off
- Chair control
- Switch buttons for public water or bottle




Arms for assistant table

Arm with adjustable height



Extended arm




Simple short arm



Optional accessories on column

Tray on simple arm

- 1x silicone tray
- 1x stainless tray
- 2x silicone tray
- 2x stainless tray




Monitor



Tray on simple arm + monitor

- 1x silicone tray
- 1x stainless tray
- 2x silicone tray
- 2x stainless tray



OP lights

FARO ALYA LED

adjustable light intensity 3,000 - 50,000 lux



FARO MAIA LED

adjustable light intensity 3,000 - 35,000 lux



FARO EDI

adjustable light intensity 8,000 - 25,000 lux



Dr LED

adjustable light intensity 3,000 lux – 40,000 lux



Foot controls

Multifunctional foot pedal



multifunctional turning



foot control (only treatment unit control)



foot switch (only instruments activation)



Patient chairs

SK 1-01



SK 1-08

