










**Since 1935**  
**CHIRANA - European quality**



**Made in Slovakia**



-  Plastic parts, electronic CHIRANA, micromotors CHIRANA, metallic parts
-  Amalgam separator, raw material, machine equipment,
-  Micromotors, scalers, upholstery for patient chair, raw material, machine equipment
-  OP lights, curing lamp, valves, syringes, suction systems
-  Driving units for patient chairs
-  Machine equipment
-  Aluminium casting parts
-  Scalers
-  Scalers

# CHIRANA CHEESE E



# Dentist tables

Dentist table 6N



Dentist table 4N



# System of instruments tubes leading

Whip system



Hanging system





# Control panels

## Button-control panel




- Spray and power set up
- Light for instrument on/off
- Chip blower function on/off
- Reverse of MM on/off
- Giromatic function on/off
- Mode control ( by foot control or manual from keyboard )
- Disinfection programm - activation
- Glass filling
- Spittoon rinsing
- Stop function
- OP light switch on/off
- Chair control

## Touch-screen control panel – only at 6N/4N dentist table available



- All functions from button-control panel
- Torque of MM adjustment with automatic stop when set torque is exceeded
- Choice of gear ratios of micromotor ( only for brushless MM )
- Program selection (only for brushless micromotor)
- Indication for instrument treatment with oil Smiloil

# Additional tray under dentist table with side control panel



1x silicone tray  
1x silicone tray with handle



1x stainless tray  
1x stainless tray with handle

2x silicone tray  
2x silicone tray with handle



2x stainless tray  
2x stainless tray with handle



# Spittoon block

## Suction systems:

- Ejector saliva
- Ejector suction
- CATTANI
- DURR MSBV
- DURR CS1
- DURR CAS1

## Hygiene:

- Clean water bottle
- Disinfection system
- Clean water bottle + disinfection system



# Arms for assistant table

Standard assistant arm 330 mm



Extended arm



Arm with adjustable height





# Assistant tables


EASY 2N – two positions



EASY 4N – four positions



## Control panel on assistant table



- Glass filling function
- Spittoon rinsing
- Stop function
- OP light switch on/off
- Chair control
- Switch buttons for public water or bottle



# Optional accessories on column

Tray table



Monitor



Tray table and monitor




# Arms of trays on column

Pantographic arm



Simple arm



# Trays on column

## **Silicone tray**

1x silicone tray  
1x silicone with handle



## **Stainless tray**


1x stainless tray  
1x stainless with handle



2x silicone tray  
2x silicone tray with handle



2x stainless tray  
2x stainless tray with handle



1x or 2x silicone tray with cup holder  
1x or 2x silicone tray with cup holder + flacons holder



1x or 2x stainless tray with cup holder  
1x or 2x stainless tray with cup holders + flacons holder



# OP lights

## FARO ALYA LED

adjustable light intensity 3,000 - 50,000 lux



## FARO MAIA LED

adjustable light intensity 3,000 - 35,000 lux



## FARO EDI

adjustable light intensity 8,000 - 25,000 lux



## Dr LED

adjustable light intensity 3,000 lux – 40,000 lux



# Foot controls

Multifunctional foot pedal



multifunctional turning



foot control ( only treatment unit control )



foot switch ( only instruments activation )



# Patient chairs

**SK 1-01**



**SK 1-08**

