

Since 1935 CHIRANA - European quality



Plastic parts, electronic CHIRANA, micromotors CHIRANA, metallic parts



Amalgam separator, raw material, machine equipment,



Micromotors, scalers, upholstery for patient chair, raw material, machine equipment



OP lights, curing lamp, valves, syringes, suction systems



Driving units for patient chairs



Machine equipment



Aluminium casting parts



Scalers



Scalers

Made in Slovakia



CHIRANA CHEESE



Dentist tables

Dentist table 6N



Dentist table 4N




System of instruments tubes leading

Whip system



Hanging system



Control panels

Button-control panel




- Spray and power set up
- Light for instrument on/off
- Chip blower function on/off
- Reverse of MM on/off
- Giromatic function on/off
- Mode control (by foot control or manual from keyboard)
- Disinfection programm - activation
- Glass filling
- Spittoon rinsing
- Stop function
- OP light switch on/off
- Chair control

Touch-screen control panel – only at 6N/4N dentist table available



- All functions from button-control panel
- Torque of MM adjustment with automatic stop when set torque is exceeded
- Choice of gear ratios of micromotor (only for brushless MM)
- Program selection (only for brushless micromotor)
- Indication for instrument treatment with oil Smiloil

Additional tray under dentist table with side control panel



1x silicone tray
1x silicone tray with handle



1x stainless tray
1x stainless tray with handle

2x silicone tray
2x silicone tray with handle



2x stainless tray
2x stainless tray with handle




Spittoon block

Suction systems:

- Ejector saliva
- Ejector suction
- CATTANI
- Durr MSBV
- Durr CS1
- Durr CAS1

Hygiene:

- Clean water bottle
- Disinfection system
- Clean water bottle + disinfection system



Arms for assistant table

Standard assistant arm 330 mm



Extended arm



Arm with adjustable height



Assistant tables


EASY 2N – two positions



EASY 4N – four positions



Control panel on assistant table



- Glass filling function
- Spittoon rinsing
- Stop function
- OP light switch on/off
- Chair control
- Switch buttons for public water or bottle



Optional accessories on column

Tray table



Monitor



Tray table and monitor




Arms of trays on column

Pantographic arm



Simple arm



Trays on column

Silicone tray

1x silicone tray
1x silicone with handle



Stainless tray


1x stainless tray
1x stainless with handle



2x silicone tray
2x silicone tray with handle



2x stainless tray
2x stainless tray with handle



1x or 2x silicone tray with cup holder
1x or 2x silicone tray with flacons holders + flacons holder



1x or 2x stainless tray with cup holder
1x or 2x stainless tray with flacons holders + flacons holder

OP lights

FARO ALYA LED

adjustable light intensity 3,000 - 50,000 lux



FARO MAIA LED

adjustable light intensity 3,000 - 35,000 lux



FARO EDI

adjustable light intensity 8,000 - 25,000 lux



Dr LED

adjustable light intensity 3,000 lux – 40,000 lux



Foot controls

Multifunctional foot pedal



multifunctional turning



foot control (only treatment unit control)



foot switch (only instruments activation)



Patient chair

SK 1-08

