

Since 1935 CHIRANA - European quality



- Plastic parts, electronic CHIRANA, micromotors CHIRANA, metallic parts
- Amalgam separator, raw material, machine equipment,
- Micromotors, scalers, upholstery for patient chair, raw material, machine equipment
- OP lights, curing lamp, valves, syringes, suction systems
- Driving units for patient chairs
- Machine equipment
- Aluminium casting parts
- Scalers
- Scalers

Made in Slovakia



CHIRANA CHEESE LIFT



Dentist tables


Dentist table 6N



Dentist table 5N



Dentist table 4N




Dentist tables

Dentist tables 5N on pantographic arm



Cart dentist tables 5N



System of instruments tubes leading

Whip system



Hanging system



Control panels

Button-control panel



- Spray and power set up
- Light for instrument on/off
- Chip blower function on/off
- Reverse of MM on/off
- Giromatic function on/off
- Mode control (by foot control or manual from keyboard)
- Disinfection programm - activation
- Glass filling
- Spittoon rinsing
- Stop function
- OP light switch on/off
- Chair control

Touch-screen control panel – only at 6N/4N dentist table available



- All functions from button-control panel
- Torque of MM adjustment with automatic stop when set torque is exceeded
- Choice of gear ratios of micromotor (only for brushless MM)
- Program selection (only for brushless micromotor)
- Indication for instrument treatment with oil Smiloil

Arms of dentist table

Short arm – 370 mm



Normal arm



Additional tray for dentist table

5N dentist table

fixed side tray




1x silicone tray
1x stainless tray
2x silicone tray
2x stainless tray



4N/6N dentist table

1x silicone tray
1x silicone tray with handle
2x silicone tray
2x silicone tray with handle



1x stainless tray
1x stainless tray with handle
2x stainless tray
2x stainless tray with handle



Spittoon block

Suction systems:

Ejector saliva
Ejector suction
CATTANI
DURR MSBV
DURR CS1
DURR CAS1

Hygiene:

Clean water bottle system
Disinfection system
Both systems



Without spittoon bowl



Assistant tables


2N – two places



4N – four places



Control panel on assistant table



- Glass filling function
- Spittoon rinsing
- Stop function
- OP light switch on/off
- Chair control
- Switch buttons for public water or bottle




Arms for assistant table

Arm with adjustable height



Extended arm




Simple short arm



Optional accessories for column

Tray table:



Monitor




Arms of trays on column

Pantographic arm

1 x silicone or 2x silicone tray
1x stainless or 2x stainless tray
with or w/o handle
with or w/o cup holder
with or w/o flacons holders



Simple arm



OP lights

FARO ALYA LED

adjustable light intensity 3,000 - 50,000 lux



FARO MAIA LED

adjustable light intensity 3,000 - 35,000 lux



FARO EDI

adjustable light intensity 8,000 - 25,000 lux



Dr LED

adjustable light intensity 3,000 lux – 40,000 lux



Foot controls

Multifunctional foot pedal



multifunctional turning



foot control (only treatment unit control)



foot switch (only instruments activation)



Patient chairs

SK 1-01



SK 1-08

