


EXTREME 2D 25 Liters



Packing

Imballo

Weight with box (EXTREME 2D 1,5HP) :

44 Kg - 97 Lbs

Box dimensions cm : 52,5 x 57,5 x 72 (0,21m³)

Quantity per (110x110) PALLET : 8 pc.

 Dental unit	Pneumatic		
	Electric-Pneumatic		
	Electric		
EXTREME 0,75	1	X	X
EXTREME 1,50	2	1	1
EXTREME 2,00	2	2	1

Compressor Model	Tank	Dryer	Cylinders	Air Intake	Air delivery	Pressure	Power	Phase	Noise	Weight	Dimensions
Modello Compressore	Serb.	Essicat.	Cilindri	Aria Aspirata	Aria resa	Pressione	Motore	Fase	Rumore	Peso	Dimensioni
Type	Liters	Complete	N°	l/min. - C.F.M.	l/min. - C.F.M.	Bar - Psi	HP - kW	Tip	dB	Kg - Lbs	cm (LxWxH)
EXTREME 2V 25L 0,75	25	X	2	200 - 7,0	60 - 2,1 (5 BAR) 42 - 1,5 (8 BAR)	8 - 114	0,75 - 0,5	Single	72	30,5 - 67	46 x 43 x 83
EXTREME 2D 25L 0,75	25	✓	2	200 - 7,0	60 - 2,1 (5 BAR) 42 - 1,5 (8 BAR)	8 - 114	0,75 - 0,5	Single	72	38,5 - 85	58,5 x 43 x 83
EXTREME 2V 25L 1,50	25	X	2	240 - 8,5	125 - 4,4 (5 BAR) 100 - 3,5 (8 BAR)	8 - 114	1,50 - 1,1	Single	72	33 - 73	46 x 43 x 83
EXTREME 2D 25L 1,50	25	✓	2	240 - 8,5	125 - 4,4 (5 BAR) 100 - 3,5 (8 BAR)	8 - 114	1,50 - 1,1	Single	72	41 - 90,5	58,5 x 43 x 83
EXTREME 2V 25L 2,00	25	X	2	350 - 12,4	156 - 5,5 (5 BAR) 127 - 4,5 (8 BAR)	8 - 114	2,00 - 1,5	Single	73	35,5 - 78	46 x 43 x 83
EXTREME 2D 25L 2,00	25	✓	2	350 - 12,4	156 - 5,5 (5 BAR) 127 - 4,5 (8 BAR)	8 - 114	2,00 - 1,5	Single	73	44 - 97	58,5 x 43 x 83

With automatic thermal protection

Con protezione termica automatica

MADE IN ITALY BY NARDI COMPRESSORI

Page 05



Packing

Imballo

Weight with box (EXTREME 2D 1,5HP) :

51 Kg - 112,4 Lbs

Box dimensions cm : 60 x 45 x 86 (0,23 m³)

Quantity per (110x110) PALLET : 8 pc.

 Dental unit	Pneumatic		
	Electric-Pneumatic		
	Electric		
EXTREME 0,75 CV	1	X	X
EXTREME 1,50 CV	2	1	1
EXTREME 2,00 CV	2	2	1
EXTREME 2,50 CV	4	3	2
EXTREME 3,00 CV	6	4	3

Compressor Model	Tank	Dryer	Cylinders	Air Intake	Air delivery	Pressure	Power	Phase	Noise	Weight	Dimensions
Modello Compressore	Serb.	Essicat.	Cilindri	Aria Aspirata	Aria resa	Pressione	Motore	Fase	Rumore	Peso	Dimensioni
Type	Liters	Complete	N°	l/min. - C.F.M.	l/min. - C.F.M.	Bar - Psi	HP - kW	Tip	dB	Kg - Lbs	cm (LxWxH)
EXTREME 2V 50L 0,75	50	X	2	200 - 7,0	60 - 2,1 (5 BAR) 42 - 1,5 (8 BAR)	8 - 114	0,75 - 0,5	Single	72	36,5 - 81	46 x 43 x 83
EXTREME 2D 50L 0,75	50	✓	2	200 - 7,0	60 - 2,1 (5 BAR) 42 - 1,5 (8 BAR)	8 - 114	0,75 - 0,5	Single	72	45 - 99	58,5 x 43 x 83
EXTREME 2V 50L 1,50	50	X	2	240 - 8,5	125 - 4,4 (5 BAR) 100 - 3,5 (8 BAR)	8 - 114	1,50 - 1,1	Single	72	39,5 - 87	46 x 43 x 83
EXTREME 2D 50L 1,50	50	✓	2	240 - 8,5	125 - 4,4 (5 BAR) 100 - 3,5 (8 BAR)	8 - 114	1,50 - 1,1	Single	72	48 - 106	58,5 x 43 x 83
EXTREME 2V 50L 2,00	50	X	2	350 - 12,4	156 - 5,5 (5 BAR) 127 - 4,5 (8 BAR)	8 - 114	2,00 - 1,5	Single	73	42,5 - 94	46 x 43 x 83
EXTREME 2D 50L 2,00	50	✓	2	350 - 12,4	156 - 5,5 (5 BAR) 127 - 4,5 (8 BAR)	8 - 114	2,00 - 1,5	Single	73	51 - 112	58,5 x 43 x 83
EXTREME 2V 50L 2,50	50	X	4	480 - 17,0	240 - 8,5 (5 BAR) 190 - 6,7 (8 BAR)	8 - 114	2,50 - 1,84	Single	73	52 - 115	64 x 45 x 89
EXTREME 2D 50L 2,50	50	✓	4	480 - 17,0	240 - 8,5 (5 BAR) 190 - 6,7 (8 BAR)	8 - 114	2,50 - 1,84	Single	73	60 - 132	64 x 45 x 89
EXTREME 2V 50L 3,00	50	X	4	600 - 21,2	373 - 13 (5 BAR) 330 - 12 (8 BAR)	10 - 145	3,00 - 2,2	Three	73	54 - 119	64 x 45 x 89
EXTREME 2D 50L 3,00	50	✓	4	600 - 21,2	373 - 13 (5 BAR) 330 - 12 (8 BAR)	10 - 145	3,00 - 2,2	Three + N	73	62 - 136	64 x 45 x 89

With automatic thermal protection

Con protezione termica automatica

EXTREME 2D 100 Liters



Packing

Imballo

Weight with box (EXTREME 2D 2,0HP) :

Kg 70 - 154 Lbs

Box dimensions cm : 60 x 45 x 142 (0,39 m³)

Quantity per (110x110) PALLET : 4 pc.

 Dental unit	Pneumatic		
	Electric-Pneumatic		
	Electric		
EXTREME 2,00	2	2	1
EXTREME 2,50	5	3	2
EXTREME 3,00	6	4	3

Compressor Model	Tank	Dryer	Cylinders	Air Intake	Air delivery	Pressure	Power	Phase	Noise	Weight	Dimensions
Modello Compressore	Serb.	Essicat.	Cilindri	Aria Aspirata	Aria resa	Pressione	Motore	Fase	Rumore	Peso	Dimensioni
Type	Liters	Complete	N°	l/min. - C.F.M.	l/min. - C.F.M.	Bar - Psi	HP - kW	Type	dB	Kg - Lbs	cm (LxWxH)
EXTREME 2V 100L 2,00	100	✗	2	350 - 12,4	156 - 5,5 (5 BAR) 127 - 4,5 (8 BAR)	8 - 114	2,00 - 1,5	Single	73	57 - 127	46 x 43 x 140
EXTREME 2D 100L 2,00	100	✓	2	350 - 12,4	156 - 5,5 (5 BAR) 127 - 4,5 (8 BAR)	8 - 114	2,00 - 1,5	Single	73	66 - 145	58,5 x 43 x 140
EXTREME 2V 100L 2,50	100	✗	4	480 - 17,0	240 - 8,5 (5 BAR) 190 - 6,7 (8 BAR)	8 - 114	2,50 - 1,84	Single	73	67 - 147	64 x 45 x 140
EXTREME 2D 100L 2,50	100	✓	4	480 - 17,0	240 - 8,5 (5 BAR) 190 - 6,7 (8 BAR)	8 - 114	2,50 - 1,84	Single	73	75 - 165	64 x 45 x 140
EXTREME 2V 100L 3,00	100	✗	4	600 - 21,2	373 - 13 (5 BAR) 330 - 12 (8 BAR)	10- 145	3,00 - 2,2	Three	73	69 - 152	64 x 45 x 140
EXTREME 2D 100L 3,00	100	✓	4	600 - 21,2	373 - 13 (5 BAR) 330 - 12 (8 BAR)	10- 145	3,00 - 2,2	Three + N	73	77 - 169	64 x 45 x 140

With automatic thermal protection

Con protezione termica automatica

MADE IN ITALY BY NARDI COMPRESSORI

Page 07



Packing

Imballo

Weight with box (EXTREME 3D 1,5HP) :

50 Kg - 112 Lbs

Box dimensions cm : 88 x 48 x 72 (0,30 m³)

Quantity per (110x110) PALLET : 8 pc.

 Dental unit	Pneumatic		
	Electric-Pneumatic		
	Electric		
EXTREME 0,75	1	X	X
EXTREME 1,50	2	1	1
EXTREME 2,00	2	2	1

Compressor Model	Tank	Dryer	Cylinders	Air Intake	Air delivery	Pressure	Power	Phase	Noise	Weight	Dimensions
Modello Compressore	Serb.	Essicat.	Cilindri	Aria Aspirata	Aria resa	Pressione	Motore	Fase	Rumore	Peso	Dimensioni
Type	Liters	Complete	N°	l/min. - C.F.M.	l/min. - C.F.M.	Bar - Psi	HP - kW	Tip	dB	Kg - Lbs	cm (LxWxH)
EXTREME 3V 50L 0,75	50	X	2	200 - 7,0	60 - 2,1 (5 BAR) 42 - 1,5 (8 BAR)	8 - 114	0,75 - 0,5	Single	72	37 - 82	84 x 31 x 68
EXTREME 3D 50L 0,75	50	✓	2	200 - 7,0	60 - 2,1 (5 BAR) 42 - 1,5 (8 BAR)	8 - 114	0,75 - 0,5	Single	72	44,5 - 98	84 x 44 x 68
EXTREME 3V 50L 1,50	50	X	2	240 - 8,5	125 - 4,4 (5 BAR) 100 - 3,5 (8 BAR)	8 - 114	1,50 - 1,1	Single	72	40 - 88	84 x 31 x 68
EXTREME 3D 50L 1,50	50	✓	2	240 - 8,5	125 - 4,4 (5 BAR) 100 - 3,5 (8 BAR)	8 - 114	1,50 - 1,1	Single	72	47,5 - 105	84 x 44 x 68
EXTREME 3V 50L 2,00	50	X	2	350 - 12,4	156 - 5,5 (5 BAR) 127 - 4,5 (8 BAR)	8 - 114	2,00 - 1,5	Single	73	43 - 95	84 x 31 x 68
EXTREME 3D 50L 2,00	50	✓	2	350 - 12,4	156 - 5,5 (5 BAR) 127 - 4,5 (8 BAR)	8 - 114	2,00 - 1,5	Single	73	50,5 - 112	84 x 44 x 68

With automatic thermal protection

Con protezione termica automatica