

# Patient Chair

## SK 1.01



 **CHIRANA**  
Medical STARÁ  
TURA



put your trust  
in tradition

 Smile

# Patient Chair SK 1.01

Reliable, ergonomically designed to give the dentist the best working position and comfort for the patient.

- Anatomical design of the backrest and seat section to avoid “stripping” effect
- Four programmable working positions
- Multi-articulating headrest
  - Thin profile backrest, available in three variations:
    - 1) Classic, 2) Wide, 3) Classic with lower cut out for headrest
  - Three variations of armrest:
    - 1) Left-Fixed, 2) Right – tilting forward, 3) Right – tilting backward.
- All Variations have the option to be upholstered
  - Safety switches on backrest and cantilever
  - High Quality seamless, stretch upholstery (air and steam permeable, hard wearing)
- Detachable hygienic foot cover
  - Cleaning of upholstery with antibacterial and antimycotic agents (do not use alcohol based cleaning agents)
  - Power supply positioned under the front cover of the patient chair
  - Foot control allows independent use of the SK1.01 chair and in conjunction with any other Smile unit
    - Max. lifting weight: 180 kg
    - Available in a large range of colours

- Min. height travel: 42 cm

- Trendelenberg movement available as optional extra

