

# Patient Chair

## SK 1.08



**CHIRANA**  
Medical STARÁ  
TURA

put your trust  
in tradition



Smile

# Patient Chair SK 1.08

New generation of patient chair. Maximum comfort for both dentist and patient.

- Anatomical design of the backrest and seat section to avoid "stripping" effect
- 2 sets of four programmable working positions with visual control on display
  - Electronic or mechanical control of multi-articulating headrest
  - Electronic control of backrest and seat movements
    - Three variations of armrest:
      - 1) fixed left, 2) tilting right,
      - 3) electronic controlled right
    - Safety microswitches with "return" movement to release a trapped object
      - New seamless upholstery technology – easy disinfection and cleaning
      - Cleaning of upholstery with antibacterial and antimycotic agents (do not use alcohol based cleaning agents)
      - High Quality stretch upholstery (air and steam permeable, hard wearing)
  - Power supply positioned under the front cover of the patient chair
  - Foot control allows independent use of the SK1.08 chair and in conjunction with any other Smile unit
    - Max. lifting weight: 185 kg
    - Available in a large range of colours
    - Trendelenberg position included as standard



Smile