

smile

FD



 **CHIRANA**
Medical STARÁ
TURÁ

put your trust
in tradition



Smile dental unit intended for dental and orthodontic treatment rooms.

This unit meets the most demanding requirements ergonomics for dentist work in lying position of the patient. Technical solution enables variety of modifications.

- Swivel motion of dentist and assistant table arms (left to right and back)
- Control of the unit functions via foil keyboard on dentist table
- Regulation of micromotor speed from 60 rpm
- Electronic "Giromatic" function of the micromotor
- Both light and acoustic warning for instrument treatment (lubrication)
- Setting the instrument power level via foil keyboard or proportionally by the foot control
- Chip Air activated by the foot control
- Automatic activation of Chip Air after completing work with the instrument
- Swivel porcelain spittoon bowl (towards the patient)
- Control of the chair is from three points: from the dentist table, from the assistant table and from the foot control
- "Chair Services" are under the front cover of the patient chair
- Available in a large range of colours
- For attaching to SK1.01 chair

- Assistant table with foil keyboard, big and small suction hoses and the option of an assistant syringe and curing light

- Dentist table with foil keyboard and the option of up to 5 instruments

