



# CURING LIGHT iLed

Cure 2mm resins in one second

## EXPERT IN LED CURING LIGHT

13 years professional technology, 39 patents, 7 Invention Patents  
(Official marks under 882, 300-2484)

light intensity up to **2300 mw/cm<sup>2</sup>**



Scan and login website  
for more information



WWW.MEDISERVIS.SK



### CONSTANT POWER OUTPUT

Best for insurance

### 360° ROTATABLE HEAD

Extend the possibility

### 1 SEC FOR CURING

Go quicker than ever

### ERGONOMICS DESIGN

Approach the hand precisely



- ✓ The use of disposable isolation sets, time-saving disinfection

### CONSTANT LIGHT



⚠ Teeth have shorter curing light without constant light intensity will cause incomplete solidification in the bottom of resin when power drop to about 10%.

- ✓ Constant light output. Solidification effect will not be affected by the consumption of remaining power.
- ✗ For most curing lights on the market, light intensity obviously decreases with the consumption of remaining power



Batteries produced in various conditions of quality level with various types of insulation allow various applications (e.g. mobile phone and some medical equipment) with no safety and some battery high temperature risk

### Battery

- ✓ Woodpecker lithium battery with CE certification, more reliable no explosion
- ✗ Batteries on the market without sound charge management, easily cause explosion under high temperature, unguaranteed quality.

# 1 SEC FOR CURING

Go quicker than ever

### BASIC DATA:

Battery	Li-Ion
Power output	AC 200V-240V

WWW.MEDISERVIS.SK